

# Pulsed NMR system

This is a work in progress...

## Interesting links

Some interesting links on pulsed nuclear magnetic resonance:

- [Pulsed NMR at UW](#)
- [Pulsed NMR at MSU](#)
- [The Basics of NMR](#) by Joseph P. Hornak

## Short description


The system consists of one in-phase/quadrature (I/Q) digital down-converter (DDC) and of one pulse generator.

The tunable frequency range covers from 0 Hz to 60 MHz.

The I/Q data rate is configurable and six settings are available: 25, 50, 125, 250, 500, 1250 kSPS.

## Hardware

The basic blocks of the system are shown on the following diagram:




The `projects/pulsed_nmr` directory contains three Tcl files: `block_design.tcl`, `rx.tcl`, `tx.tcl`. The code in these files instantiates, configures and interconnects all the needed IP cores.

## Software

The `projects/pulsed_nmr/server` directory contains the source code of the TCP server (`pulsed-nmr.c`) that receives control commands and transmits the I/Q data streams (up to  $4 \times 32 \text{ bit} \times 1250 \text{ kSPS} = 152 \text{ Mbit/s}$ ) to the control program running on a remote PC.

The `projects/pulsed_nmr/client` directory contains the source code of the control program (`pulsed_nmr.py`).



## Getting started with GNU/Linux

- Download **SD card image zip file** (more details about the SD card image can be found at [this link](#)).
- Copy the contents of the SD card image zip file to a micro SD card.
- Optionally, to start the application automatically at boot time, copy its `start.sh` file from `apps/pulsed_nmr` to the topmost directory on the SD card.
- Install the micro SD card in the Red Pitaya board and connect the power.
- Install required Python libraries:

```
sudo apt-get install python3-numpy python3-matplotlib python3
```

- Clone the source code repository:

```
git clone https://github.com/pavel-demin/red-pitaya-notes
```

- Run the control program:

```
cd red-pitaya-notes/projects/pulsed_nmr/client
python3 pulsed_nmr.py
```

## Building from source

The installation of the development machine is described at [this link](#).

The structure of the source code and of the development chain is described at [this link](#).

Setting up the Vitis and Vivado environment:

```
source /opt/Xilinx/Vitis/2020.2/settings64.sh
```

Cloning the source code repository:

```
git clone https://github.com/pavel-demin/red-pitaya-notes
cd red-pitaya-notes
```

Building `pulsed_nmr.bit`:

```
make NAME=pulsed_nmr bit
```

Building SD card image zip file:

```
source helpers/build-all.sh
```